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**For collection, transportation and storage of oral (nasal) swab DNA samples at room temperature**  
**室温保存**

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## 一般注意事项 (General Precautions)

- This preservation solution is for research use only. All due care and attention should be exercised in the handling of the solution
- Please follow the manufacture recommended sample collection procedure for maximum DNA sample protection and personal safety.

## 准备工作 (Preparation Before Starting)

- Make sure no food particle or beverage in the mouth. Rinse mouth with pure water and wait for 30 minutes before collecting sample
- Label the storage tube with name and date.

## 口腔（鼻腔）拭子 DNA 收集及保存步骤 (Oral (Nasal) Swab Collecting and Preservation Procedure)

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The saliva DNA sample stored in the preservation solution is compatible with most common DNA isolation methods. Charm Biotech's Just-a-Tube Swab Isolation kits provide easy and efficient procedure for maximum high-quality DNA isolation. The quality of purified DNA can be examined by gel electrophoresis. Quantity of purified DNA may be determined by UV spectrophotometer or fluorescent DNA assay. The expected yield is 0.5 – 5  $\mu$ g of DNA from a single swab. The actual yield of DNA will varied depending on individual donor status. Purified DNA is compatible with downstream DNA analysis, such as PCR, qPCR, NGS, southern blot, microarray studies.

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The methods for DNA purification from oral (nasal) swab stored in the Charm Biotech Oral (Nasal) Swab DNA preservation solution:

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## Troubleshooting

Problem	Cause	Solution
Low yield of product	Amount of collected cells was less than recommended. Poor quality of starting material	Be sure to follow the protocol for oral (nasal) swab collection and processing correctly.
Partial degradation	Sample is not mixed immediately and thoroughly with DNA Preservation Solution. Some dead oral (nasal) epithelial cells are already existed before collection	Be sure to put swab into preservation solution in short time (less than 5 minutes) and mixed immediately by inverting or shaking the storage tube more than 10 times. Use clean water to rinse oral or nasal cavities before collection to remove dead cells as much as possible