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- This preservation solution is for research use only. All due care and attention should be exercised in the handling of the solution
- Please follow the manufacture recommended sample collection procedure for maximum DNA sample protection and personal safety.

## 准备工作 (Preparation before Starting)

- Make sure no food particle or beverage in the mouth. Rinse mouth with pure water and wait for 30 minutes before collecting sample
- Label the storage tube with name and date.

## 唾液 DNA 收集及保存步骤 (Saliva DNA Collecting and Preservation Procedure)

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The saliva DNA sample stored in the preservation solution is compatible with most common DNA isolation methods. Charm Biotech's Just-a-Tube Saliva DNA Isolation kits provide easy and efficient procedure for maximum high-quality DNA isolation. The expected yield is 3 – 60 µg of DNA from one ml saliva. The actual yield of DNA will varied depending on individual donor status. Purified DNA is compatible with downstream DNA analysis, such as PCR, qPCR, NGS, southern blot, microarray studies.

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The methods for gDNA purification from saliva stored in the Wuhan Charm Biotech Saliva Preservation Solution:

- (1) **Just-a-Tube™ Genomic DNA Purification from Charm-Protect™ Saliva DNA Preservation Solution Kit**  
Kit features: Using Charm Biotech developed Solid Surface Reverse Binding(SSRB) technology, Single-microtube purification method. Easy and fast to perform with high rate of recovery of gDNA, Only 10 minutes of hands-on time for the purification, High yield of gDNA. Need centrifugation. (Please check [www.wuhancharmbiotech.com](http://www.wuhancharmbiotech.com) for the detail)
- (2) **Charm-Easy™ Genomic DNA Purification from Charm-Protect™ Saliva DNA Preservation Solution Kit**  
Kit features: Using Charm Biotech developed Solid Surface Reverse Binding(SSRB) technology, Single-microtube or 96-well plate purification method. NO need of centrifugation, vacuum, or magnetic force. The easiest procedure on the market. (Please check [www.wuhancharmbiotech.com](http://www.wuhancharmbiotech.com) for the detail)
- (3) **Charm-Pure™ gDNA Purification from Charm-Protect Saliva DNA Preservation Solution Kit**  
Kit features: Conventional Spin column method for easy gDNA purification.

## Troubleshooting

Problem	Cause	Solution
Low yield of product	Amount of saliva used was less than recommended. Poor quality of starting material	Be sure the protocol for saliva collection and processing is followed correctly.
Partial degradation	Saliva sample is not mixed immediately and thoroughly with Saliva Preservation Solution.	Be sure to spit equal volume of saliva to the collection vial and mixed thoroughly by inverting or shaking the vial more than 10 times.