## Charm-Pure™ Large-Scale Saliva Genomic DNA Purification

Store at room temperature, store Proteinase K and RNase A at 4 °C

### Product Contents and Storage

The components included with the Charm-Pure™ Large-Scale Saliva gDNA Purification are listed below. Upon receipt store all components at room temperature with the exception of Proteinase K and RNase A, which need to be stored at 4 °C. Sufficient reagents are included to perform 25 or 125 genomic DNA extractions.

Product Catalogue	JM--279-M	JM-279-L
Purification Scale	25	125
Binding Tube (BT4)	25 tubes	125 tubes
Lysis Buffer (PKLB)	11.0 ml	55.0 ml
Binding Buffer (KB6)	22.0 ml	110.0 ml
Spin-Column with Collection Tube(CC1)	50 sets	250 sets
Binding Buffer (GB5)	22.0 ml	110.0 ml
Washing Buffer I (WB1)	5.5 ml	27.5 ml
Washing Buffer II (WB2)	11.0 ml	55.0 ml
Elution Buffer (EB2)	50.0 ml	250.0 ml
Proteinase K (PK2) (20 mg/ml)	2.2 ml	11.0 ml
RNase A (RA2) (20 mg/ml)	2.2 ml	11.0 ml

### Product Description

The Charm-Pure™ Large-Scale Saliva Genomic DNA Purification kit is designed and optimized for easy and fast isolation of high-quality gDNA from 4 ml fresh saliva, frozen saliva or saliva sample stored in the Charm-Protect Saliva Collection and Storage Vial or preservation solutions from other companies. Based on Charm Biotech developed Solid Surface Reversible Binding (SSRB) technology and special silica membrane technology, 10 ML tubes coated with proprietary turbo-binders and spin-column with modified silica membranes are used to selectively capture and efficiently bind genomic DNA from saliva cell lysates. High-yield and high-quality DNA can be easily obtained from 4 ML saliva mixtures. The purified gDNA can be eluted in Elution Buffer or water. The purified ready-to-use gDNA is ideal for downstream application such as restriction digestion, primer extension, SNP analysis, PCR-free sequencing, whole genome amplification (WGA), DNA database/databank, and DNA methylation analysis.

### Additional Materials Needed

- 100 % ethanol , Isopropanol
- 1.5 ml or 2.0 ml microtubes, 15 ml Falcon conic centrifuge tube
- Bench-top centrifuge capable of spinning 1.5 – 2.0 ml microtubes, Vortex or Shaker, Incubator (air or water) or hot blocker
- Swinging-bucket centrifuge capable of spinning 15 ml tube, such as Eppendorf Centrifuge 5804, Beckman GS-6, or Sovall 6600

### General Precautions

- This kit is for research use only. All due care and attention should be exercised in the handling of the kits.
- Wear a laboratory coat, disposable gloves, and eye protection when handling reagents and tubes. Avoid ingestion and inhalation of reagents. Avoid skin contact with reagents in the kit. In case of contact, wash thoroughly with water. See Material Safety Data Sheets (MSDS) for emergency procedures in case of accidental contact or ingestion. All MSDS information is available upon request.
- Always use proper aseptic techniques to avoid nuclease contamination when working with nucleic acids and use only sterile, new pipette tips to prevent cross contamination.

## Preparation before Starting

- Prepare fresh working Binding Buffer prior to performing gDNA isolation from saliva lysates based on the number of samples processed. To make working Binding Buffer (KB6), mix 100 µl of Binding Buffer (KB6) with 400 µl 100 % Isopropanol and mix well. Prepare a master Binding Buffer solution based on (1) the number of samples processed and (2) any anticipated loss, generally 10 %, during dispensing. For working Binding Buffer (KB6), dispense 100 µl of Binding Buffer containing Isopropanol per 100 µl cell lysates (1:1 ratio). Discard the unused Binding Buffer at the end of the day.
- For JM-279-M kit, add 22 ml of 100 % ethanol to Washing Buffer I (WB1) and mix well. For JM-279-L kit, add 110 ml of 100 % ethanol to Washing Buffer I (WB1) and mix well. (Mark bottle that ethanol has been added!) – Store at room temperature and use WB1 containing ethanol within six (6) months.
- For JM-279-M kit, add 44 ml of 100 % ethanol to Washing Buffer II (WB2) and mix well. For JM-279-L kit, add 220 ml of 100 % ethanol to Washing Buffer II (WB2) and mix well. (Mark bottle that ethanol has been added!) – Store at room temperature and use WB2 containing ethanol within six (6) months.

## Experimental Procedure

### Preparing Lysates

The protocols listed below are for purification of gDNA from 4 ml of saliva preservation mixture. Adjust volume accordingly for actual volume of saliva mixture.

**Note: If the saliva samples are stored in the Charm-Protect Saliva Preservation Solution, please follow the procedure below to prepare the saliva lysates. (If you need to isolate gDNA from saliva samples without using Charm-Protect Saliva Preservation Solution, please check “[Experimental Procedure for Non-Charm Protected saliva samples](#)”).**

1. Add 80 µl RNase A (RA2) into a 15 ml (Falcon) conic centrifuge tube.
2. Mix the saliva sample mixture stored in the Charm-Protect Saliva Preservation Solution in the Saliva Collection Vial by inverting and shaking 7-8 times. Let the vial or tube sit for 3-5 minutes.
3. Transfer 4 ml saliva preservation mixture from the Saliva Collection Vial into the 15 ml tube preloaded with RNase A (RA2) and mix well by inverting the tube 7 – 8 times after closing the cap. (Note: Don't collect the particulate settlement at the bottom of the collection vial when aspirate saliva mixture.)
4. Add 80 µl Proteinase K (PK2) into the tube and mix well by inverting the tube 15 – 18 times after closing the cap..
5. Incubate the conic tube at 62 °C in an air (or water) incubator or a hot-block for 45 minutes with occasional mixing by tapping or inverting the tube. (Optional: Incubate the sample overnight at 52 °C if prefer)
6. Centrifuge the conic tube at  $\geq 2,250 \times g$  for 10 minutes to remove any particulate materials, such as food particles, undigested dead cells and fibers debris.

### Experimental Procedure for Non-Charm Protected saliva samples

If your saliva samples are fresh saliva, frozen saliva or saliva mixed with saliva preservation solution from other companies, please prepare saliva lysates as follows:

1. Add 80 µl RNase A (RA2) into a 15 ml (Falcon) conic centrifuge tube.
2. Add 400 µl Lysis Buffer (PKLB) into the tube (the volume ratio of PKLB vs. saliva sample is 1 to 9) .
3. Mix the saliva sample by inverting the saliva container 7-8 times and then let the container sit for 3-5 minutes.
4. Transfer 3.6 ml saliva or saliva mixture into the conic tube preloaded with RNase A (RA2) and PKLB. Mix well by inverting the 15 ml tube 9 – 10 times after closing the cap. (Note: Don't collect the particulate settlement at the bottom of the container when aspirate saliva mixture.)
5. Add 80 µl Proteinase K (PK2) into the tube and mix well by inverting the tube 15 – 18 times after closing the cap.
6. Incubate the tube at 62 °C in an air (water) incubator or a hot-blocker for 60 minutes with occasional mixing by tapping or inverting the tube. (Optional: Incubate the sample overnight at 52 °C if prefer)
7. Centrifuge the 15ml tube at  $\geq 2,250 \times g$  for 10 minutes to remove any particulate materials, such as food particles, undigested dead cells and fibers debris.

### Binding Products

1. Transfer 3.9 ml clarified saliva lysate into a Binding Tube (BT4) without collecting the debris pellet at the bottom of the 15 ml tube (Note: To avoid the debris pellet, you may leave about 150-260 µl solution in the conic tube when doing aspiration).
2. Add 3.9 ml working Binding Buffer (KB6) containing Isopropanol to the Binding Tube (BT4). Mix well by inverting the Binding Tube 15 – 18 times after closing the cap to obtain a homogenous solution. (Note: It is crucial to mix working Binding Buffer with saliva lysate thoroughly.)
3. Centrifuge the Binding Tube (BT4) at  $\geq 2,250 \times g$  at room temperature for 15 minutes to bind the DNA.
4. Open the cap, remove the solution from the Binding Tube by flipping the Binding Tube over a waste container and shaking briskly. Let the tube stand upside down on a piece of clean absorbent paper for 3 – 8 minutes to drain as much liquid as possible.

- Add 1.2 ml Elution Buffer (EB2) to the Binding Tube. Close the cap and vortex the tube for 1 minute and incubate the tube at 62 °C for 15 - 20 minute to dissolve the DNA pellet thoroughly. (You may vortex again after 62 °C incubation to help dissolve the pellet)
- Add 720 µl Binding Buffer (GB5) to the Binding Tube and mix well by inverting the tube 8-9 times after closing the cap.
- Add 720 µl ethanol to the mixture above and mix well by inverting the tube 12- 15 times (The total volume is about 2.64 ml)
- Transfer two 660 µl solution mixtures from the step above to two Spin-Columns with Collection Tubes (CC1) separately. Centrifuge the Spin-Columns at 8000 rpm (RCF = 5800 g) for one minute to allow all liquid pass through the membrane at room temperature.
- Discard the filtrate in the Collection Tube and reinsert the Spin-Column back into the Collection Tube.
- Repeat step 8 and 9 above once to allow all solution mixture pass through the membrane.

#### Washing Products

- Add 500 µl of Wash Buffer I (WB1) with ethanol to each Column.
- Centrifuge the Spin-Column at maximum speed ( $\geq 13000$  rpm or 16000 X g) for 1 minute at room temperature. Discard the filtrate.
- Add 500 µl of Wash Buffer II (WB2) with ethanol to each Column.
- Centrifuge the Spin-Column at maximum speed ( $\geq 13000$  rpm or 16000 X g) for 1 minute at room temperature. Discard the filtrate.
- Repeat step 3 and 4 above once for a total of two washes with Wash Buffer II (WB2).
- Centrifuge the spin-column at maximum speed for one minute to remove any residual wash buffer with ethanol. Discard collection tube.
- Insert the Spin-Column into a new clean micro-tube (self-provide). (Note: When transferring the spin-column, be careful of preventing the residual wash buffer contacting the tip of the column)

#### Eluting Genomic DNA

- Warm up the Elution Buffer (EB2) in a 1.5 ml micro-tube with a 62 °C hot block or water-bath.
- Add 150 µl Elution Buffer (EB2) onto the center of the Spin-Column, incubate the Spin-Column with collection micro-tube at 62 °C in an air incubator or a hot-blocker for at least 5 minutes.
- Centrifuge the column at maximum speed for 1 minute at room temperature. The micro-tube contains the filtrate which contains purified gDNA. (Note: You may perform a second elution with another 100 µl of Elution Buffer (EB2) if desired. This can increase the total gDNA yield by 15-30 %, but the final concentration of isolated DNA in the eluted sample is reduced because of the increased volume).
- Combine the two eluted gDNAs from the same saliva sample. The eluted gDNA may be used immediately in downstream applications. Alternatively, the eluted gDNA may be stored in the micro-tube at 4 °C for short-term storage or -20 °C for long-term storage.

#### Electrophoresis and Downstream Application

Purified gDNA can be examined by agarose gel electrophoresis. Yield can be measured with a spectrophotometer or fluorescent DNA assays. The yield of purified gDNA varied in the range of 20 – 120 µg from different saliva samples,  $A_{260/280} \geq 1.7$  (the majority is between 1.8 – 2.0). Purified gDNA can be used immediately for downstream application, such as restriction enzyme digestion, primer extension, SNP analysis, PCR, STR analysis, DNA sequencing, whole genome amplification (WGA) and other molecular manipulations.

#### Troubleshooting

Problem	Cause	Solution
Low DNA yield	Poor quality of starting material	Be sure the protocol for sample collection and preservation recommended by manuals is followed correctly.
	Incomplete lysis	If incomplete lysis is observed, extend the incubation time at 62 °C for an addition of 15 – 30 minutes.
	Variable DNA amount from different individual saliva samples	Increase the starting saliva sample volume for DNA isolation. Collect oral swap samples or more saliva samples
No PCR product	Missing Component in the PCR mixture	Be sure to add all components. Check positive control and negative control for PCR reaction.

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Product Catalogue	JM-279-M	JM-279-L
Purification Scale	25 次	125 次
Binding Tube (BT4)	25 个	125 个
Lysis Buffer (PKLB)	11.0 ml	55.0 ml
Binding Buffer (KB6)	22.0 ml	110.0 ml
Spin-Column with Collection Tube(CC1)	50 套	250 套
Binding Buffer (GB5)	22.0 ml	110.0 ml
Washing Buffer I (WB1)	5.5 ml	27.5 ml
Washing Buffer II (WB2)	11.0 ml	55.0 ml
Elution Buffer (EB2)	50.0 ml	250.0 ml
Proteinase K (PK2) (20 mg/ml)	2.2 ml	11.0 ml
RNase A (RA2) (20 mg/ml)	2.2 ml	11.0 ml

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