## Just-a-Tube™ Genomic DNA Purification from Saliva

**Store at room temperature, store Proteinase K and RNase A at 4 °C**

### Product Contents and Storage

The components included with the Just-a-Tube gDNA Purification from saliva are listed below. Upon receipt store all components at room temperature with the exception Proteinase K and RNase A, which need to be stored at 4 °C. Sufficient reagents are included to perform 50 or 250 genomic DNA extractions.

Product Cat. #	QT-275-M	QT-275-L
Purification Scale	50	250
Binding Tube (BT3)	50	250
Lysis Buffer (PKLB)	5.5 ml	27.5 ml
Binding Buffer (KB6)	10.0 ml	50.0 ml
Washing Buffer I (WB1)	5.5 ml	27.5 ml
Washing Buffer II (WB2)	18.0 ml	90.0 ml
Elution Buffer (EB2)	5.5 ml	27.5 ml
Proteinase K (PK2) (20 mg/ml)	1.1 ml	5.5 ml
RNase A (RA2) (20 mg/ml)	1.1 ml	5.5 ml

### Product Description

The Just-a-Tube™ Genomic DNA Purification from Saliva is designed and optimized for easy and fast isolation of high-quality gDNA from fresh saliva, frozen saliva or saliva sample stored in the Charm-Protect Saliva Collection and Storage Vial or preservation solutions from other companies. DNA sample purification can be processed just in a microtube. No silica membrane and magnetic beads are needed for DNA binding in the process. Based on Charm Biotech developed Solid Surface Reversible Binding (SSRB) technology, the Just-a-Tube system utilizes micro-tubes coated with proprietary turbo-binders acting to selectively capture and efficiently bind genomic DNA from saliva cell lysates. In the presence of Binding Buffer gDNA specifically interacts with the turbo binders and binds to the tubes while proteins and other contaminants will remain in the solution. Unbound material is removed in washing steps. The purified gDNA can easily be eluted in Elution Buffer or water. The purified ready-to-use gDNA is ideal for downstream application such as restriction digestion, primer extension, SNP analysis, PCR, sequencing, whole genome amplification (WGA) and DNA methylation analysis.

### Additional Materials Needed

- 100 % ethanol , Isopropanol
- 1.5 ml or 2.0 ml microtubes
- Bench-top centrifuge capable of spinning 1.5-2.0 ml microtubes
- Vortex or Shaker
- Incubator (air or water) or hot blocker

### General Precautions

- This kit is for research use only. All due care and attention should be exercised in the handling of the kits.
- Wear a laboratory coat, disposable gloves, and eye protection when handling reagents and tubes. Avoid ingestion and inhalation of reagents. Avoid skin contact with reagents in the kit. In case of contact, wash thoroughly with water. See Material Safety Data Sheets (MSDS) for emergency procedures in case of accidental contact or ingestion. All MSDS information is available upon request.
- Always use proper aseptic techniques to avoid nuclease contamination when working with nucleic acids and use only sterile, new pipette tips to prevent cross contamination.

## Preparation before Starting

- Prepare fresh working Binding Buffer prior to performing gDNA isolation from saliva lysates based on the number of samples processed. To make working Binding Buffer (KB6), mix 100 µl of Binding Buffer (KB6) with 400 µl 100 % Isopropanol and mix well. Prepare a master Binding Buffer solution based on (1) the number of samples processed and (2) any anticipated loss, generally 10 %, during dispensing. For working Binding Buffer (KB6), dispense 100 µl of Binding Buffer containing Isopropanol per 100 µl cell lysates (1:1 ratio). Discard the unused Binding Buffer at the end of the day.
- For QT-275-M kit, add 22 ml of 100 % ethanol to Washing Buffer I (WB1) and mix well. For QT-275-L kit, add 110 ml of 100 % ethanol to Washing Buffer I (WB1) and mix well. (Mark bottle that ethanol has been added!) – Store at room temperature and use WB1 containing ethanol within six (6) months.
- For QT-275-M kit, add 72 ml of 100 % ethanol to Washing Buffer II (WB2) and mix well. For QT-275-L kit, add 360 ml of 100 % ethanol to Washing Buffer II (WB2) and mix well. (Mark bottle that ethanol has been added!) – Store at room temperature and use WB2 containing ethanol within six (6) months.

## Experimental Procedure

### Preparing Lysates

Note: If the saliva samples are stored in the Charm-Protect Saliva Preservation Solution, please follow the procedure below to prepare the saliva lysates. (If you need to isolate gDNA from saliva samples without using Charm-Protect Saliva Preservation Solution, please check “Experimental Procedure for Non-Charm Protected saliva samples”).

### Experimental Procedure for Charm Protected saliva samples:

1. Add 20 µl RNase A (RA2) into a 1.5 ml or 2.0 ml microtube.
2. Mix the saliva sample mixture stored in the Charm-Protect Saliva Preservation Solution in the Saliva Collection Vial by inverting and shaking 7-8 times. Let the vial or tube sit for 3-5 minutes.
3. Transfer 1 ml saliva preservation mixture from the Saliva Collection Vial into the microtube preloaded with RNase A (RA2) and mix well by pipetting the solution 8-9 times or by inverting the microtube 7-8 times after closing the cap. (Note: Don't collect the particulate settlement at the bottom of the collection vial when aspirate saliva mixture.)
4. Add 20 µl Proteinase K (PK2) into the tube and mix well by inverting the tube 15-18 times.
5. Incubate the microtube at 62 °C in an air incubator or a hot-blocker for 45 minutes with occasional mixing by tapping or inverting the tube. (Optional: Incubate the sample overnight at 52 °C if prefer)
6. Centrifuge the microtube at maximum speed (about 13000 rpm or 16000 g) for 4 minutes to remove any particulate materials, such as food particles, undigested dead cells and fibers debris.

### Experimental Procedure for Non-Charm Protected saliva samples:

If your saliva samples are fresh saliva, frozen saliva or saliva mixed with saliva preservation solution from other companies, please prepare saliva lysates as follows:

1. Add 20 µl RNase A (RA2) into a 1.5 ml or 2.0 ml microtube.
2. Add 100 µl Lysis Buffer (PKLB) into the tube (the volume ratio of PKLB vs. saliva sample is 1 to 9) .
3. Mix the saliva sample by inverting the saliva container 7-8 times and then let the container sit for 3-5 minutes.
4. Transfer 900 µl saliva or saliva mixture into the microtube preloaded with RNase A (RA2) and PKLB, mix well by pipetting the solution up and down 8 - 9 times or by inverting the microtube 9-10 times after closing the cap. (Note: Don't collect the particulate settlement at the bottom of the container when aspirate saliva mixture.)
5. Add 20 µl Proteinase K (PK2) into the tube and mix well by inverting the tube 15-18 times.
6. Incubate the microtube at 62 °C in an air incubator or a hot-blocker for 45 minutes with occasional mixing by tapping or inverting the tube. (Optional: Incubate the sample overnight at 52 °C if prefer)
7. Centrifuge the microtube at maximum speed (about 13000 rpm or 16000 g) for 4 minutes to remove any particulate materials, such as food particles, undigested dead cells and fibers debris.

### Binding Products

1. Transfer 900 µl clarified saliva lysate into a Binding Tube (BT3) without collecting the debris pellet at the bottom of the microtube (Note: To avoid the debris pellet, you may leave about 150-250 µl solution in the microtube when doing aspiration).
2. Add 900 µl working Binding Buffer (KB6) containing Isopropanol to the Binding Tube (BT3). Mix well with the saliva lysate by pipetting up and down the solution 10 times or by inverting the microtube 15-18 times after closing the cap to obtain a homogenous solution. (Note: It is crucial to mix working Binding Buffer with saliva lysate thoroughly.)
3. Centrifuge the Binding Tube (BT3) at maximum speed (about 13000 rpm or 16000 g) at room temperature for 4 minutes to bind the DNA. (Note: During centrifugation, always position microtube hinge pointed outward from the center of rotation. Majority of DNA will collect at the bottom along the hinge side of the microtube.)

- Open the cap. Remove the liquid from the Binding Tube. Choose one of the methods listed below to remove the liquid. (1) Decant the solution by flipping the Binding Tube over a waste container and shaking briskly. (2) Remove the solution by aspirating the solution from the center of the bottom with a pipet tip. Be sure not to scrape the walls of the tube with pipet tips during aspiration as the products are bound to the walls.

#### Washing Products

- Add 500  $\mu$ l Washing Buffer I (WB1) containing ethanol to each tube. Mix by shaking the tube 4-5 times. Incubate at room temperature for 1 minute.
- Centrifuge the Binding Tube at maximum speed ( $\geq$  13000 rpm or 16000 g) at room temperature for one minute. (Note: During centrifugation, always position microtube hinge pointed outward from the center of rotation).
- Remove the Washing Buffer (WB1) from the Binding Tube (BT3) using one of the methods suggested in step 4 of "Binding Products".
- Add 800  $\mu$ l Washing Buffer II (WB2) containing ethanol to the Binding Tube (BT3). Mix by shaking the tube 4-5 times or by inverting the microtube 7-8 times after closing the cap.
- Centrifuge the Binding Tube at maximum speed ( $\geq$  13000 rpm or 16000 g) at room temperature for one minute. (Note: During centrifugation, always position micro-tube hinge pointed outward from the center of rotation).
- Remove the Washing Buffer (WB2) from the Binding Tube (BT3) using one of the methods suggested in step 4 of "Binding Products".
- Repeat steps 4 to 6 above once for a total of two washes with Washing Buffer (WB2).
- After the final wash, to ensure complete removal of Washing Buffer, you may spin the Binding Tube (BT3) very briefly after removing the majority of solution and aspirate the last drop of liquid at the bottom center of the tube with a 200  $\mu$ l pipet tip, and air-dry the Binding Tube (BT3) in a lab hood for 8-10 minutes to remove any residual liquid. Alternatively, put the tube in a 62 °C air incubator or a hot blocker and air-dry the tube without cover for 3-5 minutes.

#### Eluting Genomic DNA

Genomic DNA attached to the walls of the Binding Tube can be stored directly in the well for at least 6 months. If the gDNA is to be eluted, follow the procedure below

- Add 40-150  $\mu$ l Elution Buffer (EB2) into Binding Tube. (If higher concentration is preferred, add 40  $\mu$ l Elution Buffer (EB2) incubate at 62 °C for 5-10 minutes.
- Vortex the tube for 30-60 seconds or shake the tube for 5 minutes. If the liquid is cloudy or has any particulate materials inside, centrifuge the tube at maximum speed for one minute to transfer upper clear liquid to a new clean tube.
- The eluted gDNA may be used immediately in downstream applications. Alternatively, the eluted gDNA may be stored in the Binding Tube at 4 °C for short-term storage or -20 °C for long-term storage.

#### Electrophoresis and Downstream Application

Purified gDNA can be examined by agarose gel electrophoresis. Yield can be measured with a spectrophotometer or fluorescent DNA assays. A260/A280 values are in the range of 1.6-1.9 (which are lower than standard values, but there is no any problem for any downstream molecular assays based on our data showed). The yield is about 4-25  $\mu$ g. Purified gDNA can be used immediately for downstream application, such as restriction enzyme digestion, primer extension, SNP analysis, PCR, STR analysis, DNA sequencing, whole genome amplification (WGA) and other molecular manipulations.

#### Troubleshooting

Problem	Cause	Solution
Low DNA yield	Poor quality of starting material	Be sure the protocol for sample collection and preservation recommended by manuals is followed correctly.
	Incomplete lysis	If incomplete lysis is observed, extend the incubation time at 62 °C for an addition of 15-30 minutes.
	Variable DNA amount in different individual saliva samples	Increase the starting saliva sample volume for DNA isolation. Please call us for exact procedure if you want to process whole 4 ml saliva mixtures.
No PCR product	Missing Component in the PCR mixture	Be sure to add all components. Check positive control and negative control for PCR reaction.

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Product Cat. #	QT-275-M	QT-275-L
Purification Scale	50 次	250 次
Binding Tube (BT3)	50 个	250 个
Lysis Buffer (PKLB)	5.5 ml	27.5 ml
Binding Buffer (KB6)	10.0 ml	50.0 ml
Washing Buffer I (WB1)	5.5 ml	27.5 ml
Washing Buffer II (WB2)	18.0 ml	90.0 ml
Elution Buffer (EB2)	5.5 ml	27.5 ml
Proteinase K (PK2) (20 mg/ml)	1.1 ml	5.5 ml
RNase A (RA2) (20 mg/ml)	1.1 ml	5.5 ml

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©2016 武汉昌美生物 (Wuhan Charm Biotech.) Confidential, all rights reserved. More information from