## Charm-Pure™ Genomic DNA Purification from Saliva (Spin- Column)

**Store at room temperature, store Proteinase K and RNase A at 4 °C**

### Product Contents and Storage

The components included with the Charm Pure™ gDNA Purification from saliva are listed below. Upon receipt store all components at room temperature with the exception of Proteinase K and RNase A, which need to be stored at 4 °C. Sufficient reagents are included to perform 50 or 250 sample gDNA extractions.

Product Catalogue	CT-236-M	CT-236-L
Purification Scale	50	250
Spin-Column with Collection Tube (CC1)	50	250
Lysis Buffer (PKLB)	1.65 ml	8.25 ml
Binding Buffer (GB3)	8.5 ml	42.5 ml
Washing Buffer I (WB1)	5.5 ml	27.5 ml
Washing Buffer II (WB2)	11.0 ml	55.0 ml
Elution Buffer (EB2)	5.5 ml	27.5 ml
Proteinase K(PK2) (20 mg/ml)	0.55 ml	2.75 ml
RNase A (RA2) (20 mg/ml)	0.55 ml	2.75 ml

### Product Description

The Charm-Pure™ gDNA Purification from Saliva Kit is designed and optimized for easy and fast isolation of high-quality gDNA from fresh saliva, frozen saliva or saliva sample stored in the Charm-Protect Saliva Collection and Storage Vial or preservation solutions from other companies. Charm-Pure™ system employs modified silica membranes as selective nucleic acid captures in a Spin-Column format to efficiently extract cell gDNA from saliva with higher yield and purity than conventional silica membrane. DNA molecules can specifically interact with the modified silica membrane in the extraction buffer, while proteins and other contaminants are not bound and remain in the solution and passed through. After washing the column with washing buffer to remove unbound material, the purified DNA can be easily eluted in Elution Buffer or water. The purified ready-to-use gDNA is ideal for downstream application such as restriction digestion, primer extension, SNP analysis, PCR, sequencing, whole genome amplification (WGA) and DNA methylation analysis.

### Additional Materials Needed

- 100 % ethanol
- Isopropanol
- 1.5 ml or 2.0 ml microtubes
- Pipettors and pipet tips
- Bench-top centrifuge capable of spinning 1.5-2.0 ml microtubes
- Vortex or Shaker
- Incubator (air or water) or hot block

### General Precautions

- This kit is for research use only. All due care and attention should be exercised in the handling of the kits.
- Wear a laboratory coat, disposable gloves, and eye protection when handling reagents and tubes. Avoid ingestion and inhalation of reagents. Avoid skin contact with reagents in the kit. In case of contact, wash thoroughly with water. See Material Safety Data Sheets (MSDS) for emergency procedures in case of accidental contact or ingestion. All MSDS information is available upon request.

- Always use proper aseptic techniques to avoid nuclease contamination when working with nucleic acids and use only sterile, new pipette tips to prevent cross contamination.

### **Preparation before Starting**

- For CT-236-M kit, add 22 ml of 100 % ethanol to Washing Buffer I (WB1) and mix well. For CT-236-L kit, add 110 ml of 100 % ethanol to Washing Buffer I (WB1) and mix well (Mark bottle that ethanol has been added!) – Store at room temperature and use WB1 containing ethanol within six (6) months.
- For CT-236-M kit, add 44 ml of 100 % ethanol to Washing Buffer II (WB2) and mix well. For CT-236-L kit, add 220 ml of 100 % ethanol to Washing Buffer II (WB2) and mix well (Mark bottle that ethanol has been added!) – Store at room temperature and use WB2 containing ethanol within six (6) months.
- Please following the standard Charm-Protect Saliva DNA Preservation Solution collection and preservation procedure to minimize sample degradation before proceeding gDNA purification procedure in the current protocol.

## **Experimental Procedure**

### **Preparing Lysates**

The protocols listed below are for purification of gDNA from 300 µl of saliva preservation mixture. For purification of 600 µl of saliva preservation mixture, please call us for details.

Note: If the samples are Charm-Protect saliva, prepare lysates as follows (if you need to isolate other saliva samples, please see [“Procedure for Non-Charm-Protected Saliva”](#)).

#### **Procedure for saliva with Charm-Protect saliva preservation solution:**

1. Add 10 µl RNase A (RA2) into a 1.5 ml or 2.0 ml microtube.
2. Mix the saliva sample mixture stored in the Charm-Protect Saliva Preservation Solution in the Saliva Collection Vial by inverting and shaking 7-8 times. Let the vial sit for 3-5 minutes.
3. Transfer 300 µl saliva preservation mixtures from the Saliva Collection Vial into the microtube preloaded with RNase A (RA2) and mix well by pipetting the solution 8-9 times or by inverting the microtube 7 – 8 times after closing the cap. (Note: Don't collect the particulate settlement at the bottom of the collection vial when aspirate saliva mixture.)
4. Add 10 µl Proteinase K (PK2) into the saliva mixture and mix well by inverting the microtube 15 -18 times after closing the cap.
5. Incubate the microtube at 62 °C in an air incubator or a hot-blocker for 45 minutes with occasional mixing by tapping or inverting the tube. (Optional: Incubate the sample overnight at 52 °C if prefer)

#### **Procedure for Non-Charm-Protected saliva:**

If your saliva samples are fresh saliva, frozen saliva or saliva with saliva preservation solutions from other companies, please prepare lysates as follows:

1. Add 10 µl RNase A (RA2) into a 1.5 ml or 2.0 ml microtube.
2. Add 30 µl Lysis Buffer (PKLB) (the volume ratio of PKLB and saliva sample is 1 to 9).
3. Transfer 270 µl saliva or saliva preservation mixture into the microtube preloaded with RNase A (RA2) and PKLB above and mix well by pipetting the solution 8-9 times or by inverting the microtube 7-8 times after closing the cap.
4. Add 10 µl Proteinase K (PK2) into the saliva mixture and mix well by inverting the microtube 15 -16 times after closing the cap.
5. Incubate the microtube at 62 °C in an air incubator or a hot-blocker for 45 minutes with occasional mixing by tapping or inverting the tube. (Optional: Incubate the sample overnight at 52 °C if prefer.)

### **Binding DNA**

1. Add 150 µl Binding Buffer (GB3) into the saliva mixture, and add 200 µl isopropanol into the saliva mixture and mix well by pipetting the solution 8-9 times
2. Transfer all saliva solution mixture (about 670 µl) from step above to a Spin-Column with a Collection Tube (CC1).
3. Centrifuge the Spin-Column at 8000 rpm (RCF = 5800 g) for 1 minute at room temperature. Discard the filtrate in the Collection Tube and reinsert the Spin-Column back into the Collection Tube.

### **Washing DNA**

1. Add 500 µl of Wash Buffer I (WB1) with ethanol to the Column.
2. Centrifuge the Spin-Column at maximum speed (≥ 13000 rpm or 16000 g) for 1 minute at room temperature. Discard the filtrate.

3. Add 500 µl of Wash Buffer II (WB2) with ethanol to the Column.
4. Centrifuge the Spin-Column at maximum speed ( $\geq 13000$  rpm or 16000 g) for 1 minute at room temperature. Discard the filtrate.
5. Repeat step 3 and 4 above once for a total of two washes with Wash Buffer II (WB2).
6. Centrifuge the spin-column at maximum speed for one minute to remove any residual wash buffer with ethanol. Discard collection tube. And reinsert the Spin-Column into a new clean microcentrifuge tube (self-provide). (Note: When transferring the Spin-Column, be careful of preventing the residual wash buffer contacting the tip of the column)

### Eluting Genomic DNA

1. Warm up the Elution Buffer (EB2) in a 1.5 ml micro-tube with a 62 °C hot block or water-bath.
2. Add the hot 50 µl Elution Buffer (EB2) into the center of the Spin-Column and then incubated them at 62 °C for 3-5 minutes.
3. Centrifuge the column at maximum speed for 1 minute at room temperature. The microtube contains the filtrate which contains purified gDNA. (Note: You may perform a second elution with another 50 µl of Elution Buffer (EB2) if desired. This can increase the total gDNA yield by 15-30 %, but the final concentration of isolated DNA in the eluted sample is reduced because of the increased volume).
4. The eluted gDNA may be used immediately in downstream applications. Alternatively, the eluted gDNA may be stored in the microtube at 4 °C for short-term storage or -20 °C for long-term storage.

### Electrophoresis and Downstream Application

Purified gDNA can be examined by agarose gel electrophoresis. Yield can be measured with a spectrophotometer, fluorescent DNA assays or other quantification methods. For electrophoresis, loading 5-8 µl purified DNA is recommended. The purified DNA is suitable for use in restriction digestion, primer extension, SNP analysis, PCR, STR analysis, DNA sequencing, whole genome amplification (WGA) and other molecular manipulations.

### Troubleshooting

Problem	Cause	Solution
Low DNA yield	Poor quality of starting material	Be sure the protocol for sample collection and preservation recommended by manuals is followed correctly.
	Incomplete lysis	If incomplete lysis is observed, extend the incubation time at 62 °C for an addition of 15 – 30 minutes.
No PCR product	Missing Component in the PCR mixture	Be sure to add all components. Check positive control and negative control for PCR reaction.

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