Charm-Pure™ Medium-Scale Genomic DNA Purification from Saliva

Store at room temperature, store Proteinase K and RNase A at 4 °C

Product Contents and Storage

The components included with the Charm-Pure™ Medium-Scale gDNA Purification from Saliva are listed below. Upon receipt store all components at room temperature with the exception of Proteinase K and RNase A, which need to be stored at 4 °C. Sufficient reagents are included to perform 50 or 250 genomic DNA extractions.

Product Catalogue	QC-276-M	QC-276-L
Purification Scale	50	250
Binding Tube (BT3)	50	250
Spin-Column with Collection Tube(CC1)	50	250
Lysis Buffer (PKLB)	5.5 ml	27.5 ml
Binding Buffer (KB6)	10.0 ml	50.0 ml
Binding Buffer (GB5)	11.0 ml	55.0 ml
Washing Buffer I (WB1)	5.5 ml	27.5 ml
Washing Buffer II (WB2)	11.0 ml	55.0 ml
Elution Buffer (EB2)	24 ml	120 ml
Proteinase K (PK2) (20 mg/ml)	1.1 ml	5.5 ml
RNase A (RA2) (20 mg/ml)	1.1 ml	5.5 ml

Product Description

The Charm-Pure™ Medium-Scale Genomic DNA Purification from Saliva kit is designed and optimized for easy and fast isolation of high-quality gDNA from fresh saliva, frozen saliva or saliva sample stored in the Charm-Protect Saliva Collection and Storage Vial or preservation solutions from other companies. Based on Charm Biotech developed Solid Surface Reversible Binding (SSRB) technology and special silica membrane technology, micro-tubes coated with proprietary turbo-binders and spin-column with modified silica membranes are used to selectively capture and efficiently bind genomic DNA from saliva cell lysates. High-yield and high-quality DNA can be easily obtained from 1 ML saliva mixtures. The purified gDNA can be eluted in Elution Buffer or water. The purified ready-to-use gDNA is ideal for downstream application such as restriction digestion, primer extension, SNP analysis, PCR, sequencing, whole genome amplification (WGA) and DNA methylation analysis.

Additional Materials Needed

- 100 % ethanol , Isopropanol
- 1.5 ml or 2.0 ml microtubes
- Bench-top centrifuge capable of spinning 1.5-2.0 ml microtubes, Vortex or Shaker, Incubator (air or water) or hot blocker

General Precautions

- This kit is for research use only. All due care and attention should be exercised in the handling of the kits.
- Wear a laboratory coat, disposable gloves, and eye protection when handling reagents and tubes. Avoid ingestion and inhalation of reagents. Avoid skin contact with reagents in the kit. In case of contact, wash thoroughly with water. See Material Safety Data Sheets (MSDS) for emergency procedures in case of accidental contact or ingestion. All MSDS information is available upon request.
- Always use proper aseptic techniques to avoid nuclease contamination when working with nucleic acids and use only sterile, new pipette tips to prevent cross contamination.

Preparation before Starting

- Prepare fresh working Binding Buffer prior to performing gDNA isolation from saliva lysates based on the number of samples processed. To make working Binding Buffer (KB6), mix 100 µl of Binding Buffer (KB6) with 400 µl 100 % Isopropanol and mix well. Prepare a master Binding Buffer solution based on (1) the number of samples processed and (2) any anticipated loss, generally 10 %, during dispensing. For working Binding Buffer (KB6), dispense 100 µl of

Binding Buffer containing Isopropanol per 100 µl cell lysates (1:1 ratio). Discard the unused Binding Buffer at the end of the day.

- For QC-276-M kit, add 22 ml of 100 % ethanol to Washing Buffer I (WB1) and mix well. For QC-276-L kit, add 110 ml of 100 % ethanol to Washing Buffer I (WB1) and mix well. (Mark bottle that ethanol has been added!) – Store at room temperature and use WB1 containing ethanol within six (6) months.
- For QC-276-M kit, add 44 ml of 100 % ethanol to Washing Buffer II (WB2) and mix well. For QC-276-L kit, add 220 ml of 100 % ethanol to Washing Buffer II (WB2) and mix well. (Mark bottle that ethanol has been added!) – Store at room temperature and use WB2 containing ethanol within six (6) months.

Experimental Procedure

Preparing Lysates

The protocols listed below are for purification of gDNA from 1ml of saliva preservation mixture. Adjust volume accordingly for actual volume of saliva mixture.

Note: If the saliva samples are stored in the Charm-Protect Saliva Preservation Solution, please follow the procedure below to prepare the saliva lysates. (If you need to isolate gDNA from saliva samples without using Charm-Protect Saliva Preservation Solution, please check “[Procedure for Non-Charm Protected saliva samples](#)”).

Procedure for Non-Charm Protected saliva samples:

1. Add 20 µl RNase A (RA2) into a 1.5 ml or 2.0 ml microtube.
2. Mix the saliva sample mixture stored in the Charm-Protect Saliva Preservation Solution in the Saliva Collection Vial by inverting and shaking 7-8 times. Let the vial or tube sit for 3-5 minutes.
3. Transfer 1 ml saliva preservation mixture from the Saliva Collection Vial into the microtube preloaded with RNase A (RA2) and mix well by pipetting the solution up and down 8-9 times or by inverting the microtube 7-8 times after closing the cap. (Note: Don't collect the particulate settlement at the bottom of the collection vial when aspirate saliva mixture.)
4. Add 20 µl Proteinase K (PK2) into the tube and mix well by inverting the tube 15-18 times after closing the cap..
5. Incubate the microtube at 62 °C in an air incubator or a hot-blocker for 45 minutes with occasional mixing by tapping or inverting the tube. (Optional: Incubate the sample overnight at 52 °C if prefer)
6. Centrifuge the microtube at maximum speed (about 13000 rpm or 16000 g) for 4 minutes to remove any particulate materials, such as food particles, undigested dead cells and fibers debris.

Procedure for Non-Charm Protected saliva samples:

If your saliva samples are fresh saliva, frozen saliva or saliva mixed with saliva preservation solution from other companies, please prepare saliva lysates as follows:

1. Add 20 µl RNase A (RA2) into a 1.5 ml or 2.0 ml microtube.
2. Add 100 µl Lysis Buffer (PKLB) into the tube (the volume ratio of PKLB vs. saliva sample is 1 to 9).
3. Mix the saliva sample by inverting the saliva container 7-8 times and then let the container sit for 3-5 minutes.
4. Transfer 900 µl saliva or saliva mixture into the microtube preloaded with RNase A (RA2) and PKLB, mix well by pipetting the solution up and down 8-9 times or by inverting the microtube 9-10 times after closing the cap. (Note: Don't collect the particulate settlement at the bottom of the container when aspirate saliva mixture.)
5. Add 20 µl Proteinase K (PK2) into the tube and mix well by inverting the tube 15 -18 times.
6. Incubate the microtube at 62 °C in an air incubator or a hot-blocker for 45 minutes with occasional mixing by tapping or inverting the tube. (Optional: Incubate the sample overnight at 52 °C if prefer)
7. Centrifuge the microtube at maximum speed (about 13000 rpm or 16000 g) for 4 minutes to remove any particulate materials, such as food particles, undigested dead cells and fibers debris.

Binding Products

1. Transfer 900 µl clarified saliva lysate into a Binding Tube (BT3) without collecting the debris pellet at the bottom of the microtube (Note: To avoid the debris pellet, you may leave about 150-250 µl solution in the micro-tube when doing aspiration).
2. Add 900 µl working Binding Buffer (KB6) containing isopropanol to the Binding Tube (BT3). Mix well with the saliva lysate by pipetting up and down the solution 10 times or by inverting the microtube 15-18 times after closing the cap to obtain a homogenous solution. (Note: It is crucial to mix working Binding Buffer with saliva lysate thoroughly.)
3. Centrifuge the Binding Tube (BT3) at maximum speed (about 13000 rpm or 16000 g) at room temperature for 4 minutes to bind the DNA. (Note: During centrifugation, always position micro-tube hinge pointed outward from the center of rotation. Majority of DNA will collect at the bottom along the hinge side of the microtube.)
4. Open the cap, remove the solution from the Binding Tube by flipping the Binding Tube over a waste container and shaking briskly. Let the tube stand upside down on a piece of clean absorbent paper for 1-5 minutes to drain as much liquid as possible.

- Add 300 µl Elution Buffer (EB2) to the Binding Tube. Vortex for 1 minute and incubate the tube at 62 °C for 10-15 minute to dissolve the DNA pellet thoroughly.
- Add 180 µl Binding Buffer (GB5) to the Binding Tube and mix well by pipetting the solution 8-9 times.
- Add 180 µl ethanol to the mixture above and mix well by pipetting the solution 8-9 times.
- Transfer all solution mixture (about 660 µl) from the step above to a Spin-Column with a Collection Tube (CC1). Centrifuge the Spin-Column at 8000 rpm (RCF = 5800 g) for one minute to allow all liquid pass through the membrane at room temperature.
- Discard the filtrate in the Collection Tube and reinsert the Spin-Column back into the Collection Tube.

Washing Products

- Add 500 µl of Wash Buffer I (WB1) with ethanol to the Column.
- Centrifuge the Spin-Column at maximum speed (≥ 13000 rpm or 16000 g) for 1 minute at room temperature. Discard the filtrate.
- Add 500 µl of Wash Buffer II (WB2) with ethanol to the Column.
- Centrifuge the Spin-Column at maximum speed (≥ 13000 rpm or 16000 g) for 1 minute at room temperature. Discard the filtrate.
- Repeat step 3 and 4 above once for a total of two washes with Wash Buffer II (WB2).
- Centrifuge the Spin-Column at maximum speed for one minute to remove any residual wash buffer with ethanol. Discard collection tube.
- Insert the Spin-Column into a new clean microtube (self-provide). (Note: When transferring the Spin-Column, be careful of preventing the residual wash buffer contacting the tip of the column)

Eluting Genomic DNA

- Warm up the Elution Buffer (EB2) in a 1.5 ml microtube with a 62 °C hot block or water-bath.
- Add 75 µl Elution Buffer (EB2) onto the center of the Spin-Column, incubate the spin-column with collection microtube at 62 °C in an air incubator or a hot-blocker for 5 minutes.
- Centrifuge the column at maximum speed for 1 minute at room temperature. The microtube contains the filtrate which contains purified gDNA. (Note: You may perform a second elution with another 50 µl of Elution Buffer (EB2) if desired. This can increase the total gDNA yield by 15-30 %, but the final concentration of isolated DNA in the eluted sample is reduced because of the increased volume).
- The eluted gDNA may be used immediately in downstream applications. Alternatively, the eluted gDNA may be stored in the microtube at 4 °C for short-term storage or -20 °C for long-term storage.

Electrophoresis and Downstream Application

Purified gDNA can be examined by agarose gel electrophoresis. Yield can be measured with a spectrophotometer or fluorescent DNA assays. The yield of purified gDNA varied between 4-25 µg from different saliva samples, $A_{260/280} \geq 1.7$ (the majority is between 1.8-2.0). Purified gDNA can be used immediately for downstream application, such as restriction enzyme digestion, primer extension, SNP analysis, PCR, STR analysis, DNA sequencing, whole genome amplification (WGA) and other molecular manipulations.

Troubleshooting

Problem	Cause	Solution
Low DNA yield	Poor quality of starting material	Be sure the protocol for sample collection and preservation recommended by manuals is followed correctly.
	Incomplete lysis	If incomplete lysis is observed, extend the incubation time at 62 °C for an addition of 15 – 30 minutes.
	Variable DNA amount in different individual saliva samples	Increase the starting saliva sample volume for DNA isolation. Please call us for exact procedure if you want to process whole 4 ml saliva mixtures.
No PCR product	Missing Component in the PCR mixture	Be sure to add all components. Check positive control and negative control for PCR reaction.

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Product Catalogue	QC-276-M	QC-276-L
Purification Scale	50 次	250 次
Binding Tube (BT3)	50 个	250 个
Spin-Column with Collection Tube(CC1)	50 套	250 套
Lysis Buffer (PKLB)	5.5 ml	27.5 ml
Binding Buffer (KB6)	10.0 ml	50.0 ml
Binding Buffer (GB5)	11.0 ml	55.0 ml
Washing Buffer I (WB1)	5.5 ml	27.5 ml
Washing Buffer II (WB2)	11.0 ml	55.0 ml
Elution Buffer (EB2)	24 ml	120 ml
Proteinase K (PK2) (20 mg/ml)	1.1 ml	5.5 ml
RNase A (RA2) (20 mg/ml)	1.1 ml	5.5 ml

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